

Testing Donated Seeds

Since we cannot guarantee what our avid gardeners have given us, we recommend you do the following to test for viability:

Wet some kitchen paper towel so that it is damp.

Put 10 seeds on the paper towel and cover.

Place the seeds and paper towel in a zip lock bag on top of the fridge and note the seed name and date on the paper towel and the outside of the bag.

Check every few days and see if any seeds have sprouted

If all 10 seeds sprout within the germination period then the seeds are 100% viable.

Any seeds that have germinated can be potted up and grown on as seedlings.

